

Operating Instructions

***JETI* Spectroradiometer**

spectraval 15x1



JETI Technische Instrumente GmbH
Göschwitzer Straße 48
D-07745 Jena
Tel.: +49 3641 23292 00
Fax: +49 3641 23292 01
E-mail: sales@jeti.com
Internet: www.jeti.com

Contents

1	Introduction	7
2	Security Recommendations	8
3	Drivers under Windows 8.1/ 10	9
3.1	How to Install	9
3.2	How to Uninstall	10
3.3	How to Determine Driver Version	11
3.4	How to Open the Device Manager in Windows 8.1/ 10	11
4	Installation of Hardware	12
4.1	Installation of USB Device	12
4.2	Installation of Bluetooth Device spectraval	13
4.3	Installation of LAN Device spectraval 1501	14
5	Operation of Device	16
5.1	spectraval 15x1	16
5.2	Battery charging	17
5.3	Trigger Function	17
5.3.1	Connectors	17
5.3.2	External Control of the Measurement	18

6	Display Operation for spectraval 1511	19
6.1	Connection Through USB	19
6.2	Taking a Measurement	19
6.3	Display reset	20
6.4	Displaying of Results	21
6.4.1	Numeric Values	21
6.4.1.1	Numeric Values Display Options (List of Values)	21
6.4.2	Spectrum	22
6.4.2.1	Spectrum Display Options (Wavelength Range)	22
6.4.3	xy Diagram	23
6.4.3.1	xy Diagram Display Options (Color Space & White Point)	23
6.4.4	u'v' Diagram	24
6.4.4.1	u'v' Diagram Display Options (Color Space & White Point)	24
6.4.5	CRI	25
6.4.5.1	CRI Diagram Display Options (15th Color Sample)	25
6.4.6	RGB	26
6.4.6.1	RGB Display Options (RGB Color Space)	26
6.4.7	Transmission / Reflectance	27
6.4.7.1	Transm./Reflect. Display Options (Wavelength Range)	27
6.5	SD Card	28
6.5.1	Saving of Data and File Management	28
6.5.1.1	General information	28
6.5.1.2	Settings	28
6.5.2	File Explorer	29
6.6	Reference Measurement	29
6.7	Settings	31
6.7.1	Device Settings	31
6.7.1.1	Adjust brightness	31
6.7.1.2	Power Options	32
6.7.1.3	Time and date	32
6.7.1.4	Touchscreen Calibration	33
6.7.1.5	Rotate display	34
6.7.2	Measurement Settings	34
6.7.2.1	Average	35
6.7.2.2	Continuous Scan Interval	35
6.7.2.3	Sync Frequency	35
6.7.2.4	Change Calibration	36
6.7.3	View Settings	37
6.7.3.1	Units	38
6.7.3.2	Number Format	38
6.7.3.3	Show Reference as	39
6.7.3.4	Content Selection	39



7	Technical Data	40
7.1	Mechanical Dimensions	41
8	Declarations	43
8.1	Certificate of Warranty	43
8.2	CE – Declaration of Conformity	44
8.3	Declaration of UKCA Conformity	45
8.4	License Agreement	46
9	Service	48

1 Introduction

spectraval 1501 is an easy to use spectroradiometer for various applications. Its miniaturized shape makes them suitable for laboratory as well as industrial use. The units can be used in the field in combination with a laptop.

They measure the spectral Radiance of light sources from 380 nm to 780 nm within an field of view of 1.8°. Additionally photometric, colorimetric (CCT, chromaticity, color purity, dominant wavelength, CRI, CQS, RGB) and application specific data (TLCI, PAR, circadian) are calculated from this spectrum. Therefore they are suited for the test of monitors as LCD, OLED and plasma type, for the measurement of projectors and of video walls as well as for scientific investigations. The indication of the measuring spot on the target is done with a pilot laser circle.

spectraval 15x1 can be used via Bluetooth or alternatively via USB.

spectraval 15x1 can also be equipped with an optional diffusor for spectral Irradiance measurements.

Acquisition and analysis of data is done by the PC-software *JETI LiVal*. It is also possible to create customer specific software using SDK (Software Development Kit: a package of spectrometric and radiometric DLLs) provided with the instrument. It is included on the supplied USB flash drive with detailed helping information. Furthermore, the instruments can be driven by their SCPI compatible firmware commands. The command list is also available on the USB flash drive.

PC Requirements:

- Processor (x86/x64) with 1GHz or higher, min. 2 GB RAM
- Graphic resolution: 1024×768, 32 bit color depth recommended
- Windows 8.1/10

Scope of Delivery:

- Basic device with protection cap
- Operating instructions
- Tripod
- Transport box
- USB flash drive with PC software *JETI LiVal*, DLLs, operating instructions and firmware command list
- USB cable
- Bluetooth dongle
- USB battery charger

2 Security Recommendations

Please read the following safety precautions before any operation to avoid a possible bodily injury and prevent the device from damage. In order to avoid any contingent danger, this device has only to be used within the ranges specified.

- The target pointer is a class 1 laser product¹, which has to be handled carefully. Being directed into human eyes, it can cause loss of eyesight. Do not look into the laser beam or direct the laser to another person!
- Avoid using liquid cleaners.
- Never spill liquids of any kind on the device.
- Never touch the glass window of the device, nor diffusers, glasses or filters contained in additional equipment.
- If any mechanical damages on the device are detected, please contact *JETI*.

The device must be disposed of in accordance with the applicable regulations. It contains one (spectraval 1501) or two (spectraval 1511) lithium ion batteries.

¹Classification according to IEC 60825-1

3 Drivers under Windows 8.1/ 10

3.1 How to Install

Do **not** connect the device during driver installation!

Note: You must have administrator privileges to install the device driver.

- Insert the USB flash drive delivered with your device in an open USB slot on your computer.
- Run install.exe from the USB flash drive. The following window will appear:



- Click the **Install USB drivers...** button and follow the instructions to install the drivers.
- If a warning window appears that 'Windows can't verify the publisher of this software', click on **Install Anyway**.
- After finishing, connect your device to an USB-port. Windows completes the driver installation automatically.


If the driver could not be installed successfully, please uninstall the driver according to instructions from the next chapter, and try again.

3.2 How to Uninstall

This chapter describes how to uninstall old driver software.

Note: You must have administrator privileges to uninstall the device driver.

- Run 'CDMUninstallerGUI.exe'. You can find this program on your *JETI* USB flash drive in the directory Tools ► CDMUninstaller.



- Enter '0403' in the Vendor ID text box. Leave the Product ID text box blank. Click on **Add** to add the list entry.
- Repeat the previous step and type '0C6C' in the Vendor ID text box and add it to the list in the same way.
- Make sure that no *JETI* USB-device is connected to the computer and click **Remove Devices**.
- Please wait while the drivers are uninstalling. This may take some time.
- The system should now be cleaned from old drivers. Close the program by clicking on **Cancel**.

3.3 How to Determine Driver Version

Note: For this procedure, the device, which was already installed, must be connected to the computer.

JETI USB driver consists of two parts, a device driver and a virtual COM-port driver. For the device to function correctly, the versions of the two drivers need to be the same.

To determine the version of the virtual COM-port driver:

- Open the 'device manager' (see next chapter on how to do this).
- Open the branch 'Ports (COM & LPT)' by clicking on the small plus-sign.
- Right-click the **JETI specbos xxxx Port** (or USB Serial Port) entry and choose **Properties** from the context menu.
- In the appearing property window open the tab **Driver** by clicking on it.
- On the appearing page you can find the driver version and the driver's date.

To determine the version of the device driver:

- Open the 'device manager' (see next chapter on how to do this).
- Open the branch 'Universal Serial Bus controllers' by clicking the small plus-sign.
- Right-click the **JETI specbos xxxx Device** (or USB Serial Converter) entry and choose **Properties** from the context menu.
- In the appearing property window open the tab **Driver** by clicking on it.
- On the appearing page you can find the driver version and release date.

3.4 How to Open the Device Manager in Windows 8.1/ 10

The Windows *Device Manager* shows information about the devices installed on your computer.

There are several ways to open the *Device Manager*, depending on the Windows version and the settings of Windows. The easiest of them is:

- Press **Win + Break** on the keyboard.
- Click **Device Manager** on the left side of the appearing window.


4 Installation of Hardware

4.1 Installation of USB Device


Unpack the spectroradiometer **spectraval 15x1** carefully and check the delivered parts.

Attention: The USB driver must be installed at first without connected instrument.

The **spectraval** devices are powered by USB if operated via USB interface. Additionally, the battery will be charged during this time.



Control elements of **spectraval 1501**



Control elements of **spectraval 1511**

Connect the device to a USB interface of a PC or laptop using the enclosed USB cable. Shortly afterwards the device is ready for operation, e.g. by *JETI LiVal*.

Reset: Please take a little rod, e. g. of a paper chip and press the reset button through the hole in the device housing.

4.2 Installation of Bluetooth Device spectraval

spectraval can be used via Bluetooth.


The Bluetooth interface is often integrated into modern computers. If your computer is not equipped with it, you need to use the Bluetooth dongle included in the delivery or any other one. Please connect it to the computer and switch on the **spectraval** device.

The software guides you through the installation.


- Click on the Bluetooth symbol.




- Click Add a Bluetooth Device.



- After a short time the spectraval device will be recognized. Click Ready to pair.



- If you are asked for a pin number, use 0000.
- If the following screen appears, the installation is complete and the device is ready to use:



- Close the window.

4.3 Installation of LAN Device spectraval 1501



The **spectraval** device and the controlling PC have to use the same network address space, to communicate with the **spectraval-LAN** via a network connection.


Remark: The **spectraval 1501-LAN** needs a 5 V power supply if it is used via network or USB.

The default network setting of the **spectraval** device is set to DHCP (obtain network settings automatically). If a DHCP server is running in your network, the device will obtain an IP address from within the local network address space, and no further steps are required.

If no DHCP server is available (e.g. during a direct connection to a PC) or you want to give the device a fixed IP address for some other reasons, the network setting of the **spectraval** device has to be changed through a tool called 'JETI_LAN_settings.exe' (available on the installation USB flash drive).

Note: You must have administrator privileges!

- Connect **spectraval 1501-LAN** via USB to computer.
- Start 'JETI_LAN_settings.exe' from Tools ▶ JETI_LAN_Settings to change LAN configurations.
- A window will appear where you can configure the network setting of the **spectraval** device.



- After clicking on **Set** button a message appears that device parameters have been changed successfully.

Remark: The IP address of a PC can be found using 'cmd.exe' in the command line and then entering the command 'ipconfig'.


5 Operation of Device

5.1 spectraval 15x1


The instrument will be switched on automatically shortly after it was connected to a PC when the USB cable is used. This will be shown by the shining green LED.

Attention: The USB driver must be installed at first without connected instrument (see chapter 3 'Drivers under Windows 8.1/ 10', page 9).

If you communicate via Bluetooth, the device has to be switched on using the 'On/ off' switch on the rear side.



spectraval 1501



spectraval 1511

Remark: Connect the device to your PC to charge the battery via USB.

If **spectraval** is connected to PC via Bluetooth and the connection is interrupted, **spectraval** is switched off immediately.

The battery status is shown by the yellow LED on the rear side of the device.

Status LEDs

Green:	Device is on
Yellow:	On Device is charging Off Device is not charging Flashing Charging error
Red:	Battery low. Device will be switched off shortly after this.
Blue:	Active Bluetooth connection

5.2 Battery charging

The battery can be charged using the USB battery charger or using the USB port of a PC. To use USB port for charging, plug in the USB device cable and make sure that the device is recognized by Windows correctly.

The battery will remain in the instrument, it is not exchangeable by the user.




If the device is operated by *JETI LiVal* software, battery status can be seen in the title bar of the window and through menu **Info** > **Battery Status**. If battery is low, corresponding warnings will be shown.

5.3 Trigger Function

5.3.1 Connectors

It is possible to trigger the measurement externally with a trigger female connector at the rear side of the device. It has the following pin out:



A suitable male connector is included in the device's delivery package.

5.3.2 External Control of the Measurement

A measurement can be initiated externally by a shortcut between the ground (GND) and trigger (TRG) pins. This can be done, e.g., by a foot push-button. Furthermore, it is possible to use a TTL signal (3 V to 5 V) to start a measurement. This is done with the falling edge of the signal.


Both triggering versions act similarly as pressing the green Measurement button of the *JETI LiVal* software.

See the firmware description to change the settings.

6 Display Operation for spectraval 1511

6.1 Connection Through USB

Once a **spectraval 1511** is connected via USB, the following window is displayed:




To use the **spectraval 1511** as a stand alone device through its display functionality, tap on the 'Play' button at the lower-left corner. To use it with *JETI LiVal*, tap on Trigger or press the hardware Measurement button on the device: this will cause *JETI LiVal* to initiate a measurement and receive data from the device.

6.2 Taking a Measurement

This section describes how to get started by taking a quick measurement. Note that some of the references to buttons and lights are describes in the previous sections.

- Turn on the **spectraval 1511** by pressing the 'On/Off' switch. The green status LED will illuminate.
- The **spectraval 1511** will initialize and show the main measuring interface:



- The measuring area can now be marked using the **Target** button or or by briefly pressing the 'On/ off' switch.
- Aim the laser point by pointing the measuring head on the front of the device in the direction of light. Press the **Measurement** button.
- The results will be instantly displayed on the screen.


6.3 Display reset

With **spectraval 1511** it can happen that the display crashes and nothing works anymore. Then press the 'On/Off' switch and immediately after the measuring button. Wait a few seconds and calibrate the display.

6.4 Displaying of Results


6.4.1 Numeric Values

This screen shows a list of measurements recorded after pressing the Measurement button as demonstrated in the previous section.



6.4.1.1 Numeric Values Display Options (List of Values)


The Numeric Values list can be customized with different units of measure (e.g. Lv, xy, CCT) through the context-based 'Display options', accessed from the Numeric Values Screen.



Tap on the checkbox to select/deselect an item from the list and then return to the Numeric Values by tapping the Back icon. Now your new data-set is displayed.


6.4.2 Spectrum

Light Spectrum of the latest measurement taken.



6.4.2.1 Spectrum Display Options (Wavelength Range)


The start and the end of the displayed spectrum can be adjusted through the context-based 'Display options', accessed from the Spectrum Screen.



This setting is common with the setting for the transmission/reflectance diagram.


6.4.3 xy Diagram

The CIE1931 chromaticity x, y values of the latest measurement taken.



6.4.3.1 xy Diagram Display Options (Color Space & White Point)


By tapping the checkbox you can display the RGB color triangle and the white point of the chosen RGB color Space, through the context-based 'Display options'.



These settings are common with the settings of the u'v' diagram.


6.4.4 u'v' Diagram

The CIE1976 chromaticity u' , v' values of the latest measurement taken.



6.4.4.1 u'v' Diagram Display Options (Color Space & White Point)


By tapping the check box you can display the RGB color triangle and the white point of the chosen RGB color Space, through the context-based 'Display options'.



These settings are common with the settings of the xy diagram.


6.4.5 CRI

The CRI R_1 – R_{15} values of the latest measurement taken.




6.4.5.1 CRI Diagram Display Options (15th Color Sample)

By tapping the checkbox you can select whether the 15th CRI color sample ('Asian skin color') is used, through the context-based 'Display options'.




6.4.6 RGB

The RGB values of the latest measurement taken.




6.4.6.1 RGB Display Options (RGB Color Space)

The required RGB color space can be selected (e.g. PAL, NTSC, HDTV) through the context-based 'Display options', accessed from the RGB Screen.




6.4.7 Transmission / Reflectance

Transmission / Reflectance curve. A reference must be measured or chosen (see chapter 6.6 'Reference Measurement', page 29).



6.4.7.1 Transm./Reflect. Display Options (Wavelength Range)

The start and the end of the displayed spectral data can be adjusted through the context-based 'Display options', accessed from the Transm./Reflect. Screen. This is the same wavelength range which is used by Spectrum Screen.



This setting is common with the setting for the spectrum diagram.

6.5 SD Card

spectraval 1511 can store measuring data on SD cards (only support 1 GB up to 32 GB Fat32).

6.5.1 Saving of Data and File Management

6.5.1.1 General information

spectraval 1511 supports two file formats: '.csv' for all usual measurements and for reference files and '.tr' for transmission/reflectance data.

Every measurement can be saved by tapping on the **Save** button from the main window. If you are in Transmission/Reflectance mode, a '.tr' file will be produced. In all other cases a '.csv' file will be produced.

The **Save** button in the Transmission/Reflectance mode is specially marked to avoid ambiguity:




6.5.1.2 Settings

Go to **File Settings** in the Main Settings window (alternatively long tap on the **Save** icon in the main window).

Use the long button beneath the 'Automatically Save Each Measurements' checkbox to select a directory your measurements should be saved to.

Names are given to files automatically according to some pattern. There are two options: if the checkbox 'Use Time Stamp as File Prefix' is set, time stamp is used as file names; if this checkbox is not set, some user specified (under **Edit File Prefix**) prefix with a number is used.


'Automatic Save' feature allows to write every measurement with a name specified by a certain pattern as described above onto the SD card automatically. If automatic saving is on, **Save** button in the main window is marked red.



6.5.2 File Explorer

spectraval 1511 has a simple file explorer which allows to perform some basic file operations: cut, copy, paste, rename, create a new directory, delete.


To find it: go to Main Settings / File Settings (or long-tap on **Save** button in the main window) and then tap on **File Explorer**.



6.6 Reference Measurement


The Reference item allows you to compare two measurement side by side: the current measurement against historical data from SD card.

If you selected or measured reference in any way described below, reference data will appear (see also chapter 6.7.3.3 'Show Reference as', page 39) and Transmission/Reflectance curve will become available.






Generally there are two options:

- If you want to use some file you already have saved on the SD card as reference, long-tap on the REF-field in the top line and select a '.csv' file. Confirm your choice with the **Select** button.



- If you don't intend a long-term usage of the references and want to renew it before every measuring session, short tap on the REF-field in the top line. 'Press 'Meas.' text appears. If you do a measurement now, it will be saved automatically to a hidden temporary file and it will be set as reference.






6.7 Settings

For a better understanding of the possible settings some important terms and measurement parameters are explained below. The different adjustment possibilities are divided into **Device Settings**, **Measurement Settings** and **View Settings**.

6.7.1 Device Settings




6.7.1.1 Adjust brightness

The brightness of the display can be changed under option Adjust brightness. Just click on the bar (full control corresponds to 100 % brightness)

6.7.1.2 Power Options

The Power Setting Options allow you to set different actions to take when the device is left idle for a certain amount of time.




Do not forget to click **Save** button.

If the device will turn off after a specific time, any setting will be lost!

6.7.1.3 Time and date



When the device is delivered, the CET is set. In this menu you can adjust the time.



Do not forget to click **Save** button.

6.7.1.4 Touchscreen Calibration



Touchscreen Calibration opens a window in which you can recalibrate the touchscreen (realign) should it react imprecisely to touch. At first you are asked, if you really want to calibrate.



What you need to do is to touch consequently 3 points that will be displayed.


6.7.1.5 Rotate display

Rotate display allows a 180° rotation of the display. At first you are asked, if you really want to rotate.



What you need to do is to touch consequently 3 points that will be displayed. Please calibrate the display after changed the alignment.

6.7.2 Measurement Settings



6.7.2.1 Average


The Average setting can be used to obtain the spectrum as the mean value from up to 254 successive measurements.

6.7.2.2 Continuous Scan Interval

If you want to measure continuously, you can select the interval in the range of 0 s to 65 535 s. Note that you can begin a measurement session by tapping the Measurement button in the different Mode View and stop the session by pressing the button again.


6.7.2.3 Sync Frequency

If you want to measure modulated light sources (typical examples are TV monitors and PWM modulated LED's), it is necessary to synchronize your device to the repetition rate, because **spectraval** devices are very sensitive. Otherwise the results, especially of bright sources, will fluctuate significantly. If the repetition rate of the display lies in the order of magnitude of the measuring time of the instrument, the repeatability of the measurement will be poor.



To do this, go to Measuring Settings via Settings-icon of the main screen, or just long-tap on the SYNC-field in the top line. There you have two options: it is possible to

set the frequency manually between 16 Hz and approximately 2000 Hz or to let the device to determine it. If latter use the button Determine.



If the sync measurement fails, it can have the following reasons:

- The instrument is positioned too far from the light source.
- The influence of the surrounding light or other light sources is too significant.
- The source is not modulated (note that some monitors do not fluctuate in intensity if their brightness control is set to maximum).


6.7.2.4 Change Calibration

Additional accessories are recognized by the device automatically via I²C interface. The software will use automatically the appropriate calibration file if the attached accessory is compatible with the device. By default, the first position is reserved for the Luminance mode measurement without an accessory.


Remark: Always make sure that the measuring mode shown in the MODE-field corresponds the attached accessory. Otherwise the device will use a wrong calibration file and the measuring results may be confusing.

JETI spectraval 15x1 are delivered in a calibrated status. The calibration data are saved in the device. A separate file exists for each measuring mode and each accessory.

The list of all calibrations saved in the device is available through the **View** button right from 'List of Calibrations' under Measuring Settings. Or just long-tap on the MODE-field in the main screen.




6.7.3 View Settings




6.7.3.1 Units

Lets to switch between SI and imperial units. In case of imperial units, Luminance units are footlamberts (*fL*) and Illuminance units are footcandle (*fc*).




6.7.3.2 Number Format

Allows the user to switch between scientific and floating point notation.




6.7.3.3 Show Reference as

Lets to choose the appearance of reference values: as absolute values or as differences between measured values and reference values.



6.7.3.4 Content Selection

If you want to make a selection that is different to the default, select/deselect a mode view from the list. Now only the selected modes will be displayed.



7 Technical Data

Spectral range	380 nm to 780 nm
Optical resolution (FWHM)	4.5 nm
Wavelength accuracy	± 0.5 nm
Calculated wavelength step	1 nm
Digital electronic resolution	16 bit ADC
Measuring values	Spectral Radiance, integral Luminance and Radiance, Chromaticity x, y; u', v'; Correlated color temperature; Dominant wavelength, color purity; Color rendering index CRI; circadian metrics; RGB; CQS; TLCl; etc
Measuring range Luminance	0.2 cd/m ² to 180 000 cd/m ² (CIE source A) typical 0.2 cd/m ² to 140 000 cd/m ² (typical warm white LED)
Viewing angle	1.8°
Measuring distance/ diameter in luminance mode	20 cm - $\varnothing$ 6 mm; 100 cm - $\varnothing$ 31 mm
Luminance accuracy	± 4.4 % (CIE source A @ 100 cd/m ² , k=2)
Luminance repeatability	± 1 % (CIE source A)
Color chromaticity accuracy	± 0.002 x, y (CIE source A, k=2)
Color reproducibility	± 0.0005 x, y (CIE source A)
CCT reproducibility	± 20 K (CIE source A)
Integration time range	10 μ s ... 60 s
Dispersive element	Diffraction grating
Light receiving element	CCD line array with 2048 pixels
Operating conditions	Temperature 10 °C to 40 °C Humidity <85 % relative humidity at 35 °C
Power supply	USB powered 5 V power supply (LAN-version only)
Pilot laser	Laser Class 1, $\lambda = 645 \dots 660$ nm (classific. accord. to EN60825-1/Oct. 2003)

PC interface

1501	USB 2.0 fullspeed (12 Mbit/s); Bluetooth optional LAN
1511	USB 2.0 fullspeed (12 Mbit/s); Bluetooth

Dimensions


spectraval 1501	140×80×70mm
spectraval 1511	140×115×70mm

Weight


spectraval 1501	400 g
spectraval 1511	500 g

7.1 Mechanical Dimensions

The following drawings show the position of the thread downside the device and the general dimensions of the instrument.



spectraval 1501



spectraval 1511

8 Declarations

8.1 Certificate of Warranty

Certificate of warranty for

Model: **Spectroradiometer**

Type: **spectraval 15x1**

JETI Technische Instrumente GmbH (referred to as *JETI*) hereby warrants this equipment as follows:

If any part of this unit (except as described below) fails due to poor workmanship or material (determined by *JETI*) within TWO (2) years from date of delivery, that part will be exchanged at no charge.

This warranty is valid only when the unit is installed and adjusted according to factory specifications and serviced by competent authorized personnel.

JETI does not assume responsibility for any of the following, all of which are excluded from the coverage of this warranty:

- damage due to ordinary wear and tear, abusive use, or lack of proper maintenance;
- damage due to harsh mechanical shock, e.g. falling to the floor;
- loss or damage due to adverse environmental conditions or acts of God;
- loss of wages or income due to repair, replacement, malfunction or damage.

Warranty becomes void if serial number is removed or defaced, or the instrument was opened by the customer.

If a defect appears which the customer feels is covered by this warranty, a written notice describing the defect must be sent to *JETI*'s office at the current address of record.

Upon receipt of customer's written report of a defect if the defective items are covered by this warranty, *JETI* will repair or replace it at no charge to the customer, within 30 days after receipt of the returned unit (provided there are no labour problems or materials shortages to cause delays). The choice between repair and replacement is to be made by *JETI* based on actual conditions or circumstances. *JETI* reserves the right to substitute new and improved equipment or parts at any time. The obligation to replace defective parts does not require replacement of the complete unit.

NO OTHER WARRANTY EXPRESSED OR IMPLIED IS APPLICABLE TO THIS UNIT.

This warranty is in effect for a period of 24 months, beginning one week after the date of delivery.

8.2 CE – Declaration of Conformity

No: **1605J1029V1.0**

We:

JETI Technische Instrumente GmbH Göschwitzer Straße 48 07745 Jena GERMANY
--

declare that the product series

JETI spectroradiometer **spectraval 15x1**

to which this declaration relates is in conformity with the requirements of following documents:

- EU Electromagnetic Compatibility (EMC) Directive 2014/30/EU

EU Harmonised standards:

- EN 61326-1 (2018-09)
- EN 55011 (2017-03)
- EN 61000-4-2 (2009-12); EN 61000-4-3 (2011-04)
- EN 61000-4-8 (2010-11)

The **CE** sign confirms the conformity of the product with the standards and directives mentioned above.

JETI
Technische Instrumente GmbH
Göschwitzer Straße 48, D-07745 Jena
Tel.: +49 (0)3641 232 92 00
Fax: +49 (0)3641 232 92 01



Jena, December 2019

8.3 Declaration of UKCA Conformity

We:

JETI Technische Instrumente GmbH Göschwitzer Straße 48 07745 Jena GERMANY
--

hereby declares that the product series

JETI spectroradiometer **spectraval 15x1**

placed on the market by the Company and its subsidiaries are compliant with Electromagnetic Compatibility Regulations 2016 (S:I: 2016/1091) and following documents.

EU Harmonised standards:

- EN 61326-1 (2018-09)
- EN 55011 (2017-03)
- EN 61000-4-2 (2009-12); EN 61000-4-3 (2011-04)
- EN 61000-4-4 (2013-04); EN 61000-4-6 (2014-08)
- EN 61010-1 (2020-03)

The conformity of the product with the standards and directives mentioned above, is confirmed by the UKCA sign.

JETI
Technische Instrumente GmbH
Göschwitzer Straße 48, D-07745 Jena
Tel.: +49 (0)3641 232 92 00
Fax: +49 (0)3641 232 92 01



Jena, December 2021

8.4 License Agreement

License Agreement for JETI Technische Instrumente GmbH Software Products
IMPORTANT – READ CAREFULLY BEFORE USING THE JETI SOFTWARE

JETI hereby grants to you a non-exclusive license to use one copy of the delivered JETI software program (the 'software') on a single computer according to the following items:

The software is in 'use' on the computer when it is loaded into temporary memory (i.e. RAM) or installed into permanent memory (e.g. hard disk, CD-ROM, or other storage device) of the computer. However, installation on a network server for the sole purpose of distribution to one or more other computer(s) is not constitute 'use' for which a separate license is required, provided you have separate license for each computer to which the software is distributed.

It is only allowed to use the software in connection with the JETI product (instrument, spectrometer or OEM electronics), which was delivered with this software. For additional JETI products you need an extra software license.

The software is owned by JETI or its suppliers and is protected by German copyright laws and international treaty provisions and all other applicable national laws. Therefore, you must treat the software like any other copyrighted material (e.g. a book or musical recording) except that if the software is not copy protected you may either

- make one copy of the software solely for backup or archival purposes, or
- transfer the software to a single hard disk provided you keep the original solely for backup or archival purposes. You may not copy the user documentation provided with the software and the connected JETI product.

You may not:

- permit other individuals to use the software except under the terms listed above;
- modify, translate, reverse engineer, decompile, disassemble (except to the extent that this restriction is expressly prohibited by law) or create derivative works based upon the software or documentation;
- copy the software or documentation (except for back-up or archival purposes);
- rent, lease, transfer, or otherwise transfer rights to the software or documentation;
- remove any proprietary notices or labels on the software or documentation.

Any usage in a forbidden way terminates your Software license immediately.

- You agree that you are authorized only to use the software and documentation in a manner that complies with all applicable laws in the jurisdictions in which you use the software and documentation, including, but not limited to, applicable restrictions concerning copyright and other intellectual property rights.
- You may not use the software in an attempt to, or in conjunction with, any device, program or service designed to circumvent technological measures employed to control access to, or the rights in, a content file or other work protected by the copyright laws of any jurisdiction.
- If you receive the first copy of the software electronically and a second copy on physical media (e.g., CD, diskette, etc.), the second copy may be used for archival purposes only and may not be transferred to or used by any other person. This license does not grant you any right to any enhancement or update.

The software may include certain external software components ('Add-Ons'), including the Lab Windows user interface and drivers. You may only call to or otherwise use such Add-Ons through the use of the JETI applications. Any direct use of Add-Ons through a non-JETI proprietary application, including a custom or user-written application is prohibited by this Agreement.

If you have purchased JETI software in connection with a JETI product, you have the right of purchasing a software update, if it is available, for a priority price.

Title, ownership, rights, and intellectual property of the software and documentation remain in JETI and/or its suppliers. Title, ownership rights and intellectual property rights in and to the content accessed through the software including the content of software media demonstration files are to be retained by the applicable content owner and may be protected by applicable copyright or other law. This license gives you no rights to such content.

For the JETI operating environment, the following terms apply:

JETI warrants that for a period of ninety (90) days from the date of acquisition the software, if operated as directed, will substantially achieve the functionality described in the documentation. JETI does not warrant however that your use of the software will be uninterrupted or that the operation of the software will be error-free or secure. JETI also warrants that the media containing the software, if provided by JETI, is free from defects in material and workmanship and will so remain for ninety (90) days from the date you acquire the software.

To the maximum extent permitted by applicable law JETI and its suppliers disclaim all other warranties either express or implied including but not limited to implied warranties of merchantability and fitness for a particular purpose with regard to the software, the accompanying written materials and any accompanying hardware.

If the media is subjected to accident abuse or improper use during the warranty period; or if you violate the terms of this License Agreement, this warranty is immediately terminated. This warranty is void if the software is used on or in conjunction with hardware or software other than the unmodified version of hardware and software with which the software was designed to be used as described in the documentation.

This limited warranty gives you specific legal rights.

Your remedies:

JETI's sole liability for a breach of this warranty is to be in JETI's sole discretion:

- to replace your defective media; or
- to advise you how to achieve substantially the same functionality with the software as described in the documentation through a procedure different from that set forth in the documentation; or
- if the above remedies are impracticable, to refund the license fee, if any, you paid for the software. Repaired, corrected or replaced software and documentation covered by this limited warranty for the period remaining under the warranty that covered the original software or if longer for thirty (30) days after the date *JETI* either shipped to you the repaired or replaced software or advised you as to how to operate the software so as to achieve the functionality described in the documentation, whichever is applicable. Only if you inform *JETI* of the problem with the software during the applicable warranty period and provide evidence of the date you acquired the software will *JETI* be obligated to honor this warranty.

Under no circumstances and under no legal theory whether in tort contract or otherwise *JETI*, its suppliers or resellers are liable to you or any other person for any indirect, special, incidental or consequential damages of any character including without limitation damages for loss of goodwill, work stoppage, computer failure or malfunction or any and all other commercial damages or losses even if *JETI* have been informed of the possibility of such damages or for any claim by any other party. Further, in no event *JETI*'s liability under any provision of this agreement exceeds the license fee paid to *JETI* for the software and documentation.

This software is intended for use only with properly licensed media, content, and content creation tools. It is your responsibility to ascertain whether any copyright, patent or other licenses are necessary and to obtain any such licenses to serve and/or create, compress or download such media and content. You agree to record, play back and download only those materials for which you have the necessary patent, copyright and other permissions, licenses, and/or clearances. You agree to hold harmless, indemnify and defend *JETI*, its officers, directors and employees, from and against any losses, damages, fines and expenses (including attorneys' fees and costs) arising out of or relating to any claims that you have

- viewed, downloaded, encoded, compressed, copied or transmitted any materials (other than materials provided by *JETI*) in connection with the software in violation of another party's rights or in violation of any law, or
- violated any terms of this License Agreement.

This License Agreement will automatically terminate if you fail to comply with any term hereof. No notice must be required from *JETI* to effect such termination.

Technical support for the software, as made available by *JETI*, is given directly from the *JETI* office in Jena, Germany.

This License Agreement constitutes the complete and exclusive agreement between *JETI*, its resellers and the end user, notwithstanding any variance with any purchase order or other written instrument submitted by you, whether formally rejected by *JETI* or not. The acceptance of any purchase order you place is expressly made conditional on your consent to the terms set forth herein. The terms and conditions contained in this License Agreement may not be modified except in a writing duly signed by you and an authorized representative of *JETI*. If any provision of this License Agreement is held to be unenforceable for any reason, such provision is to be reformed only to the extent necessary to make it enforceable, and such decision does not affect the enforceability of such provision under other circumstances, or of the remaining provisions hereof under all circumstances. This License Agreement is governed by the laws of Germany, without regard to conflicts of law provisions, and you hereby consent to the exclusive jurisdiction of the provincial and federal courts sitting in Germany. Any and all unresolved disputes relating in any way to, or arising out of, the software, your usage of the software or this License Agreement are to be submitted to arbitration in Germany; except that, to the extent that you have breached or have indicated your intention to breach this License Agreement in any manner which violates or may violate *JETI*'s intellectual property rights, or may cause continuing or irreparable harm to *JETI* (including, but not limited to, any breach that may impact *JETI*'s intellectual property rights, or a breach by reverse engineering), *JETI* may seek injunctive relief, or any other appropriate relief, in any court of competent jurisdiction.

Copyright © 2021 JETI Technische Instrumente GmbH and/or its suppliers.

JETI Technische Instrumente GmbH, Göschwitzer Straße 48, 07745 Jena, Germany. All rights reserved.

9 Service

Please contact in case of any question or technical problem:

JETI Technische Instrumente GmbH
Göschwitzer Straße 48
D-07745 Jena
GERMANY
Tel.: +49 3641 23292 00
Fax: +49 3641 23292 01
E-mail: support@jeti.com
Internet: www.jeti.com

Copyright © 2021 JETI Technische Instrumente GmbH. All rights reserved.
Software and operating instruction are delivered with respect to the License agreement and can be used only in accordance with this License agreement. The hard and software as well as the operating instruction are subject to change without notice.
JETI Technische Instrumente GmbH assumes no liability or responsibility for inaccuracies and errors in the operating instruction.
It is not allowed to copy this documentation or parts of it without previous written permission by
JETI Technische Instrumente GmbH

December 9, 2021