



JETI Technische Instrumente GmbH
Göschwitzer Strasse 48
D - 07745 Jena
Germany
Tel. : +49 3641 23292 00
Fax : +49 3641 23292 01
e-mail : sales@jeti.com
Internet : www.jeti.com

Technical Note 28

When I should select **spectraval 15x1** instead of **specbos 1211**?



spectraval 1501 is the replacement of **specbos 1201**. It covers the same spectral range (VIS) but is about ten times more sensitive. Therefore, it can be used for low light applications in the same manner as **specbos 1211**. **spectraval 1501** is intended for Radiance mode measurements only, a diffusor has to be ordered separately (ACC 032). If a stand alone solution is needed, then **spectraval 1511** can be used. It has the identical parameters as **spectraval 1501** but is equipped with display and SD card.

When I should use **spectraval** instead of **specbos 1211**? Please consider the following aspects:

- If a wavelength range beyond visible is not needed (version of **spectraval 15x1** up to 1 μm is planned)
- If no accessory than a diffusor is needed (currently only a diffusor is available for **spectraval**)
- If a stand-alone operation is needed (**spectraval 1511** is equipped with display and SC card slot)
- If a wireless connection is preferred (**spectraval** is equipped with Bluetooth interface by default)

The following table shows a comparison of the major parameters of all three devices:

Parameter	spectraval 1501/ 1511	specbos 1211-2	specbos 1201 (not available)
Wavelength range	380 ... 780 nm (350 ... 1000 nm)	350 ... 1000 nm UV version: 230 nm ... 1000 nm	380 ... 780 nm
Optical resolution (FWHM)	4.5 nm (2 nm - highres)	4.5 nm	5 nm
Measuring range luminance (Illum. A/ typical warm white LED spectrum)	0.2 ... 180 000 cd/m ² / 0.2 ... 140 000 cd/m ²	0.2 ... 150 000 cd/m ² / 0.2 ... 100 000 cd/m ²	2 ... 70 000 cd/m ² / 2 ... 50 000 cd/m ²
Measuring range illuminance (Illum. A/ typical warm white LED spectrum)	2 ... 1 250 000 lx / 1 ... 1 000 000 lx (with optional diffusor)	1 ... 1 800 000 lx / 1 ... 1 500 000 lx	20 ... 500 000 lx / 20 ... 370 000 lx
Viewing angle Radiance mode	1.8 °	1.8 °	1.8 °
Measuring time on grayscale pattern of a monitor	~ 0.4 s @ 100 cd/m ² ~ 3 s @ 5 cd/m ²	~ 0.5 s @ 100 cd/m ² ~ 4 s @ 5 cd/m ²	~ 2.5 s at 100 cd/m ² ~ 30 s at 5 cd/m ²
Interfaces	USB and Bluetooth Optional: LAN instead of Bluetooth	USB Optional: additional Bluetooth, RS232 or LAN	USB
Software	JETI LiVal	JETI LiVal	JETI LiVal
Power supply	USB Battery USB power supply	USB LAN: 9V power supply Bluetooth: battery, 9 V power supply	USB
Weight	1501: 400 g 1511: 500 g	450 g	350 g
Dimensions	140 mm x 80/ 115 mm x 70 mm	180 mm x 82 mm x 53 mm	140 mm x 58 mm x 34 mm

Please note: specbos 1201 is no longer available. The parameters are only listed in the table to help customers who own a device to compare more easily.